

Load data downloaded from SPRING to Partek

Rapo 2019 03 29

Overview of steps

1. Download data of interest from the SPRING^{1,2} interactive data explorer online (or run locally).
2. Run `spring_to_partek.py` in Terminal using downloaded data as input.
3. Upload data to Partek³

[1] Weinreb C, Wolock S, Klein AM. SPRING: a kinetic interface for visualizing high dimensional single-cell expression data. *Bioinformatics*. 2018 Apr 1;34(7):1246–8.


[2] <https://kleintools.hms.harvard.edu/tools/spring.html>

[3] <http://www.partek.com>

Download data from SPRING

Using the following dataset as an example⁴:

https://kleintools.hms.harvard.edu/tools/springViewer_1_6_dev.html?datasets/mouse_HPCs/basal_bone_marrow/full




[4] Tusi BK, Wolock SL, Weinreb C, Hwang Y, Hidalgo D, Zilionis R, et al. Population snapshots predict early haematopoietic and erythroid hierarchies. *Nature*. 2018 Mar 1;555(7694):54–60.

Download data from SPRING

Select cells to
download




Download



Untar the downloaded file, place in one directory
together with spring_to_partek.py

Run SPRING to partek

Open Terminal*




Navigate to the folder with spring_to_partek.py

```
cd ~/Desktop/demo
```

Run the python** code to prepare data for partek

```
python spring_to_partek.py "demo.155387179961"
```

The output is a folder called “for_partek”




* Instructions for Mac

** Written for python 3, specifically my python version is: Python 3.6.6 :: Anaconda custom (64-bit)


Upload data to Partek

New project > import data > import single cell data



Upload data to Partek


Select file

Select single cell files 

☐ Partek Flow Server ☒ My Computer ☐ URL

no files selected

No file-size limit



Select 3 files: barcodes.tsv, features.tsv, matrix.mtx
Press "Next"

Upload data to Partek

Partek Flow Queue Projects

Home > demo > Import single cell data > File format options

<input checked="" type="checkbox"/>	Sample name	Files	Cells
<input checked="" type="checkbox"/>	Project_demo_18	Project_demo_18	346

Annotation

Use annotation file ☐

Counts format

Raw counts ☐

Counts in log base

Report features without counts ☒

Gene deduplication

Deduplication method ☒ Mean ☐ Maximum ☐ Sum

Back Finish


Continue with default setting, except for the Raw counts box, counts from spring are not raw, they are normalized.

Press “Finish”.



Wait, check “Queue” for progress.

Upload data to Partek: Attributes from SPRING

Click on node
(red rectangle)



Under pre-analysis
tools, click on
Annotate cells



Browse > My computer >
choose file > select attributes.txt

Upload data to Partek: Attributes from SPRING


From here, further instructions by Partek:

<https://documentation.partek.com/display/FLOWDOC/Annotate+cells>


Attribute name	Terms	Import	Attribute type
Cell name	cell_1056, cell_1082, cell_1108,...	<input type="checkbox"/>	Categorical
numeric_P(Erythroid)	0.0, 0.00020563, 0.00024764, 0...	<input checked="" type="checkbox"/>	Numeric
numeric_P(Granulocyte)	0.00012293, 0.0015918, 0.0026...	<input checked="" type="checkbox"/>	Numeric
numeric_P(Lymphoid)	0.0, 0.00021641, 0.0003168, 0...	<input checked="" type="checkbox"/>	Numeric
numeric_P(Dendritic cell)	0.0, 0.0001238, 0.00013700000...	<input checked="" type="checkbox"/>	Numeric
numeric_P(Megakaryocyte)	0.0, 0.00010067, 0.00011106, 0...	<input checked="" type="checkbox"/>	Numeric
numeric_P(Monocyte)	0.0, 0.0001695, 0.01003499999...	<input checked="" type="checkbox"/>	Numeric
numeric_P(Basophil/Mast cell)	0.0, 0.00013663, 0.00014897, 0...	<input checked="" type="checkbox"/>	Numeric
numeric_Total counts	1005.0, 1017.0, 1031.0, 1035.0,...	<input checked="" type="checkbox"/>	Numeric
numeric_Uniform	0.0	<input checked="" type="checkbox"/>	Numeric
Louvain cluster	10, 11, 2, 3, 6, 8	<input checked="" type="checkbox"/>	Categorical
Library	1, 2, 3, 4	<input checked="" type="checkbox"/>	Categorical
Prep Batch	1, 4	<input checked="" type="checkbox"/>	Categorical

The code `spring_to_partek.py` prepends the prefix “numeric_” to colotracks that were numeric in SPRING. All other colotracks are categorical.

Example: make t-SNE plot with Partek and color by uploaded attributes



Example: make t-SNE plot with Partek and color by uploaded attributes



t-SNE by Partek, colotrack "Louvain cluster" from SPRING